

Guangzhou Roson Medical Devices Co., Ltd.

Catalogue-Nitrate Detector 2020







Add: Room 201, Floor 2, Building A Self-coding Workshop, No. 17 Hejing South Road, Dongsha Industrial Park, Liwan District, Guangzhou City, Guangdong Province, China, Post: 510385

Mobile: (+86)18998156065

E-mail: sale5rs@rosonmd.com

Contact Person: Jessica Lyu

Website: www.rosonmd.en.alibaba.com

Model	Product Picture	Product Description	Package Picture	Packing Details
HK-GT1		<p>Nitrate Detector (Fruit+Vegetable)</p> <p>Range of measurement: 20-10000mg/kg Testing time: 3 seconds Testing error: <22.5% Power: lithium battery USB charger: 250mA Power: 5v</p>		<p>Box size: 15*15*5cm Weight: 85g 1 set GT1 including 1 main machine 1 USB charger 1 bag 1 manual 1 box</p>
HK-GT2		<p>Nitrate Detector (Fruit+Vegetable+Meat)</p> <p>Range of measurement: 20-10000mg/kg Testing time: 3 seconds Testing error: <22.5% Power: lithium battery USB charger: 250mA Power: 5v</p>		<p>Box size: 15*15*5cm Weight: 85g 1 set GT1 including 1 main machine 1 USB charger 1 bag 1 manual 1 box</p>
HK-GT3		<p>Nitrate Detector (Fruit+Vegetable+Meat+Water)</p> <p>Range of measurement: 20-10000mg/kg Testing time: 3 seconds Testing error: <22.5% Power: lithium battery USB charger: 250mA Power: 5v</p>		<p>Box size: 15*15*5cm Weight: 85g 1 set GT1 including 1 main machine 1 USB charger 1 bag 1 manual 1 box</p>